

# Single-cell RNA-Seq





Bulk





Single Cell



# Single Cell RNA-Seq: a recent exploding technology



# Isolating Single Cells is a big challenge

	Micro-manipulation / Automated Pipetting	FACS	Microwell encapsulation	Droplet encapsulation
				
Cell Stress	Low	Moderate	Moderate	Moderate
Selection	Yes	Yes	No* / Yes**	No*
Doublet	Low	Low	Low-High	Moderate
Throughput	Low	Moderate	Moderate	High
Capture efficiency	Low	Moderate	Moderate	Low-Moderate
Academic / Commercial scRNA workflow	- CellenONE (Cellenion) <sup>†</sup> - Smart-Seq2 (42)	- MARS-Seq (39) - Smart-Seq2 (42)	- C1 (Fluidigm) - ddSeq (Biorad / Illumina) - iCell8 (Clontech) <sup>†*</sup> - Rhapsody (BD)	- InDrop (1 CellBio) - DropSeq (Dolomite-bio) - 10X (Chromium)
Example of use	Fragile rare cells	Rare cells based on phenotype or marking	Large cell numbers	Large cell numbers

	FACS		Microwell encapsulation				Droplet encapsulation		
	Smart-Seq2	MARS-Seq	C1	ddSeq	iCell8	Rhapsody	InDrop	DropSeq	10X
Singlet Capture efficiency	82%	92%	39%	2.6%	37%++	Not reported	7%	Not reported	50%
Doublet rate	Not reported	2.27%	3-30%	5.8%	1.3-4%	0.6%	4%	0.36-11.3	1.6-3%
Reference	42	39	37 FWP	PB	PB	PB	36	37	26


<sup>†</sup>Automated pipetting system

<sup>\*</sup>Preselection or enrichment can be performed prior

<sup>++</sup>Only reagents added to wells containing singlets, determined by system


FWP: Fluidigm white paper

PB: Product brochure / manual



<https://www.frontiersin.org/articles/10.3389/fimmu.2018.01553/full>


# Smart-seq(2)



Smart-seq (current version V4, the most used/well known is V2) is used to sequence the **full-length transcriptome** of individual cells

**Smaller throughput**, usually < 1000 cells

Since fewer cells are sequenced, sequencing depth tends to be higher and thus **more genes are captured**



<https://www.nature.com/articles/nprot.2014.006>

Requires laborious lab work  
Usually for detailed analysis of specific  
hard to get cell populations (eg. pollen)


# Chromium 10x Genomics

Uses microdroplet encapsulation for individual lysis and amplification

Current version is v3.1; has 3' and 5' sequencing variants (**not full-length transcripts**)

**High throughput** (5-10k cells/well) – good for organs or even whole small organisms


Tends to get smaller amount of RNA per cell / **noisier data**




<https://www.nature.com/articles/ncomms14049>

# A few extra notes

- Some applications add “spike-in” RNA
  - Used to control quality of samples, since we know the amounts of spike-in beforehand
  - Spikes often found identified as “ERCCxxx”
- The use of UMIs are becoming quite common
  - To avoid effect of amplification of small amounts of RNA (particularly for 10x genomics)
  - Each molecule is tagged with a UMI at library preparation
    - Counts with same UMI should be counted as only 1
  - Smart-Seq2 does not have it (smart-seq3 now includes)



# Several factors need to be taken into account



In the end, you use what you have at your institute